


Epsylar+

FOR DIAGNOSTIC AND INTERVENTIONAL
ACCESS TO THE PERIPHERAL VASCULATURE



Epsylar+

GUIDING SHEATH

Suitable for peripheral interventions, the Epsylar+ Guiding Sheath offers more possibilities as it can be used as an introducer sheath or a guiding catheter. It is designed for the introduction of diagnostic and interventional devices into the vasculature.

SMOOTH TRANSITION AND EASY INSERTION

FEATURES

Platinum-Iridium marker makes precise locating and measuring



Coiled stainless steel provides superior torqueability and support

Hydrophilic coating on distal end reduces vasospasm and smoothens access

PTFE inner layer allows devices to pass through smoothly

Silky transition between sheath and dilator minimizes vessel trauma


Atraumatic tip for safe operation





FEATURES

The detachable hemostatic valve pairs with a wide range of accessories, giving the physician greater operational flexibility



ORDER INFORMATION



Box: 1 Unit



Sheath length (cm)	French size					Tip configuration
	4F	5F	6F	7F	8F	
40	VSS4F01	VSS5F01	VSS6F01	VSS7F01	VSS8F01	straight
55	VSS4F02	VSS5F02	VSS6F02	VSS7F02	VSS8F02	
70	VSS4F03	VSS5F03	VSS6F03	VSS7F03	VSS8F03	
90	VSS4F04	VSS5F04	VSS6F04	VSS7F04	VSS8F04	

YOUR PATIENT IS THE FOCUS OF OUR EVERYDAY ACTIONS

Optimed, based in Ettlingen (Germany), has been engaged in research, development, manufacturing and worldwide distribution of high-end medical products for minimal invasive therapy since 1996. More than 200 employees in our headquarters, production facilities and our global sales teams ensure every day that we provide you and your patients with the highest quality medical devices.



Visit our website at www.optimed.com to discover our comprehensive portfolio and to view our latest product catalogue.



Order Code 9997-0034EN
Rev (08/2026) | © optimed

Distributed by:
optimed Medizinische Instrumente GmbH
Ferdinand-Porsche-Straße 11 | 76275 Ettlingen, Germany
Tel. +49 (0) 7243 / 76 33-0 | Fax +49 (0) 7243 / 76 33-624
info@optimed.com | www.optimed.com

CE 0123

Manufactured by:
Shanghai INT Medical Instruments Co., Ltd.
No. 925, Jinyuan Yi Road, 201803 | Shanghai, P.R.C.
Tel. +86 21 59140056

Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80 | 20537 Hamburg, Germany
Tel: +49 40 2513175