
optimed

Cemento-MP

TUNED SYSTEM FOR VERTEBROPLASTY






Cemento-MP

THE REUSABLE VERTEBROPLASTY SYSTEM

The vertebroplasty system consists of the **reusable application system Cemento-MP,** the special vertebroplasty needles, injection cylinder, and connection system.



FEATURES

Large screw applicator for a powerful and controlled injection

Special vertebroplasty needle featuring shock-resistant bayonet Luer lock fittings

Pressure-resistant injection cylinder with special reinforcement and rotating Luer lock

Reusable application system


The special vertebroplasty needle has indentations on the needle and stylet that indicate the opening of the bevel edge at the needle tip. The special bevel edge makes it easier to control the needle and the cement delivery.



Pressure resistant injection cylinder with special reinforcement and extra large, rotating Luer-lock




Relief lever for immediate pressure release



The injection cylinder also features special reinforcements.

Adjusting screw allows for a quick and non-messy exchange of injection cylinders



RAPID EXCHANGE CYLINDER MECHANISM



Reduction adapter – Standard Luer-Lock



Aspiration cannula

CEMENTO-MP
REUSABLE VERTEBROPLASTY SYSTEM



Cemento-MP	
	Order code
Reusable application system incl. flushing adapter	1382-0000



Cemento Mini Set for Cemento-MP	
	Order code
Supplementary set suitable for the Cemento-MP system	1382-0100



Special vertebroplasty needle set for Cemento-MP		
Diameter	Needle length	Order code
mm / gauge	mm	
3,2 / 10	100	1382-1010
	150	1382-1510
2,4 / 13	100	1382-1013
	150	1382-1513



Special vertebroplasty needle for Cemento-MP		
Diameter	Needle length	Order code
mm / gauge	mm	
3,2 / 10	100	1394-1010
	150	1394-1510
2,4 / 13	100	1394-1013
	150	1394-1513
1,8 / 15	100	1394-1015
	150	1394-1515

Order Code 9997-0019EN
Rev (09/2026) | ©optimed

Distributed by:
optimed Medizinische Instrumente GmbH
Ferdinand-Porsche-Straße 11 | 76275 Ettlingen Germany
Tel. +49 (0) 7243 / 76 33-0 | Fax +49 (0) 7243 / 76 33-624
info@optimed.com | www.optimed.com

Manufactured by:
PAJUNK® Medical Produkte GmbH
Karl-Hall-Straße 1
78187 Geisingen | Germany

CE 0124